

## Anti-Sortilin antibody (ECD) (STJ13100315)

STJ13100315

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Sortilin (ECD) is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Rat/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	IgG purified
<b>Dilution Range</b>	IHC, WB. A dilution of 1:1000 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
<b>Formulation</b>	Lyophilised
<b>Isotype</b>	IgG
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">20661</a>
<b>Gene Symbol</b>	<a href="#">Sort1</a>
<b>Uniprot ID</b>	<a href="#">SORT_MOUSE</a>
<b>Immunogen</b>	A synthetic peptide from the extracellular domain of rat Sortilin conjugated to blue carrier protein was used as the antigen.
<b>Immunogen Region</b>	ECD
<b>Specificity</b>	Specific for Sortilin.
<b>Immunogen Sequence</b>	