

Anti-Syntaxin16 antibody (STJ13100316)

STJ13100316

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Syntaxin16 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:500 to 1:1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	228960
Gene Symbol	Stx16
Uniprot ID	STX16_MOUSE
Immunogen	A synthetic peptide from mouse Syntaxin16 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in chicken, orangutan and pig.
Immunogen Region	
Specificity	Specific for Syntaxin16.
Immunogen Sequence	