

Anti-SUMO1 antibody (STJ13100334)

STJ13100334

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SUMO1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Rat/Mouse/Human/Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	7341
Gene Symbol	SUMO1
Uniprot ID	SUMO1_HUMAN
Immunogen	A synthetic peptide from mouse SUMO 1 conjugated to blue carrier protein was used as the antigen. The antigen is homologous in human, rat and marmoset.
Immunogen Region	
Specificity	Specific for SUMO1.
Immunogen Sequence	