

Anti-5HT1D R antibody (C-Term) (STJ13100353)

STJ13100353

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-5HT1D R (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Marmoset/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	3352
Gene Symbol	HTR1D
Uniprot ID	5HT1D_HUMAN
Immunogen	A synthetic peptide from c-terminal region of human 5HT1D receptor conjugated to blue carrier protein was used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for 5HT1D receptor.
Immunogen Sequence	