

Anti-SMURF1 antibody (Internal Region) (STJ13100358)

STJ13100358

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SMURF1 (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	57154
Gene Symbol	SMURF1
Uniprot ID	SMUF1_HUMAN
Immunogen	A synthetic peptide from internal region of mouse E3 ubiquitin protein ligase SMURF1 conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	Internal Region
Specificity	Specific for SMUF1.
Immunogen Sequence	