

Anti-Spinophilin antibody (450-500 aa) (STJ13100376)

STJ13100376

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Spinophilin (450-500 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 1000 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	84687
Gene Symbol	PPP1R9B
Uniprot ID	NEB2_HUMAN
Immunogen	A synthetic peptide from aa region 450-500 of human Spinophilin conjugated to blue carrier protein was used as the antigen. The peptide is homologous in mouse and rat.
Immunogen Region	450-500 aa
Specificity	Specific for Spinophilin.
Immunogen Sequence	