

Anti-TAS1R2 antibody (230-280 aa) (STJ13100437)

STJ13100437

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAS1R2 (230-280 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	83770
Gene Symbol	Tas1r2
Uniprot ID	TS1R2_MOUSE
Immunogen	A synthetic peptide from aa region 230-280 of mouse TAS1R2 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	230-280 aa
Specificity	Specific for TAS1R2.
Immunogen Sequence	