

Anti-VAMP1, 2, 3 antibody (25-65) (STJ13100474)

STJ13100474

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VAMP1, 2, 3 (25-65) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant proteins). A working concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	9341 6844 6843 VAMP3 VAMP2 VAMP3_HUMAN VAMP2_HUMAN VAMP1_HUMAN
Immunogen	A synthetic peptide from aa region 25-65 of human VAMP 3 conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	25-65
Specificity	Specific for VAMP1, 2, 3.
Immunogen Sequence	