

Anti-TRPV3 antibody (STJ13100492)

STJ13100492

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Guinea Pig polyclonal antibody anti-TRPV3 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Guinea Pig
Reactivity	Rat/Mouse/Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 1:500 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	246788
Gene Symbol	Trpv3
Uniprot ID	TRPV3_MOUSE
Immunogen	A synthetic peptide rat TRPV3 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in mouse.
Immunogen Region	
Specificity	Specific for TRPV3.
Immunogen Sequence	