

Anti-Synaptophysin antibody (STJ13100498)

STJ13100498

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Guinea Pig polyclonal antibody anti-Synaptophysin is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Guinea Pig
Reactivity	Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 500-10000 is recommended for WB and for IHC-P. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	20977
Gene Symbol	Syp
Uniprot ID	SYPH_MOUSE
Immunogen	A synthetic peptide from rat Synaptophysin conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for Synaptophysin.
Immunogen Sequence	