

Anti-VPS29 antibody (STJ13100508)

STJ13100508

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VPS29 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	51699
Gene Symbol	VPS29
Uniprot ID	VPS29_HUMAN
Immunogen	A synthetic peptide from mouse VPS29 conjugated to blue carrier protein was used as the antigen. The antigen is homologous in rat, human, chicken, zebrafish and xenopus.
Immunogen Region	
Specificity	Specific for VPS29.
Immunogen Sequence	