

Anti-VGluT3 antibody (STJ13100512)

STJ13100512

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VGluT3 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human/Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 500 is recommended for WB and 1:250 for IHC-P. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	246213
Gene Symbol	SLC17A8
Uniprot ID	VGLU3_HUMAN
Immunogen	A synthetic peptide from mouse VGLUT3 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for VGLUT3.
Immunogen Sequence	