

Anti-SLC6A19 antibody (350-400 aa) (STJ13100522)

STJ13100522

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SLC6A19 (350-400 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Marmoset

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	664630
Gene Symbol	Slc6a19
Uniprot ID	S6A19_RAT
Immunogen	A synthetic peptide from aa region 350-400 of mouse SLC6A19 conjugated to blue carrier protein was used as the antigen. The antigen is homologous in rat.
Immunogen Region	350-400 aa
Specificity	Specific for SLC6A19.
Immunogen Sequence	