

Anti-ACE2 antibody (750-800 aa) (STJ13100523)

STJ13100523

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Guinea Pig polyclonal antibody anti-ACE2 (750-800 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Guinea Pig
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 1: 500 to 1:1000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	302668
Gene Symbol	Ace2
Uniprot ID	ACE2_RAT
Immunogen	A synthetic peptide from aa region 750-800 of mouse ACE2 conjugated to blue carrier protein was used as the antigen. The antigen is homologous in rat.
Immunogen Region	750-800 aa
Specificity	Specific for ACE2.
Immunogen Sequence	