

Anti-ATG10 antibody (STJ13100527)

STJ13100527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-ATG10 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. Use at a dilution of 1:200 to 1: 2000. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	83734
Gene Symbol	ATG10
Uniprot ID	ATG10_HUMAN
Immunogen	A synthetic peptide from human ATG10 conjugated to an immunogenic carrier protein has been used as the antigen.
Immunogen Region	
Specificity	Specific for ATG10.
Immunogen Sequence	