

## Anti-ATG3 antibody (STJ13100530)

STJ13100530

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Sheep polyclonal antibody anti-ATG3 is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC/WB
<b>Host/Source</b>	Sheep
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	IHC, WB. Use at a dilution of 1:200 to 1:2000. The optimal dilution should be determined by the end user.
<b>Formulation</b>	Lyophilised
<b>Isotype</b>	
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">64422</a>
<b>Gene Symbol</b>	<a href="#">ATG3</a>
<b>Uniprot ID</b>	<a href="#">ATG3_HUMAN</a>
<b>Immunogen</b>	A synthetic peptide from human autophagy-related protein 3 (ATG3, APG3, APG3L, APG3-like, hApg3, protein PC3-96) conjugated to an immunogenic carrier protein has been used as the antigen.
<b>Immunogen Region</b>	
<b>Specificity</b>	Specific for ATG3.
<b>Immunogen Sequence</b>	