

Anti-ALMS1 antibody (4100-4150 aa) (STJ13100616)

STJ13100616

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ALMS1 (4100-4150 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	7840
Gene Symbol	ALMS1
Uniprot ID	ALMS1_HUMAN
Immunogen	A synthetic peptide from aa region 4100-4150 of human ALMS1 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	4100-4150 aa
Specificity	Specific for ALMS1.
Immunogen Sequence	