

Anti-ABCA4 antibody (1450-1500 aa) (STJ13100654)

STJ13100654

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ABCA4 (1450-1500 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	24
Gene Symbol	ABCA4
Uniprot ID	ABCA4_HUMAN
Immunogen	A synthetic peptide from aa region 1450-1500 of human ABCA4 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	1450-1500 aa
Specificity	Specific for ABCA4.
Immunogen Sequence	