

Anti-BACE1 antibody (350-400 aa) (STJ13100754)

STJ13100754

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-BACE1 (350-400 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. Use at a dilution of 1:200 to 1:2000. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	23621
Gene Symbol	BACE1
Uniprot ID	BACE1_HUMAN
Immunogen	A synthetic peptide from aa region 350-400 of human BACE1 conjugated to an immunogenic carrier protein has been used as the antigen.
Immunogen Region	350-400 aa
Specificity	Specific for BACE1.
Immunogen Sequence	