

Anti-Pan-RhO antibody (50-100 aa) (STJ13100830)

STJ13100830

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-Pan-RhO (50-100 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	6010
Gene Symbol	RHO
Uniprot ID	OPSD_HUMAN
Immunogen	A synthetic peptide from aa region of 50-100 of human RhOA conjugated to blue carrier protein was used as the antigen.
Immunogen Region	50-100 aa
Specificity	Specific for RhO family.
Immunogen Sequence	