

Anti-LL37 antibody (STJ13100856)

STJ13100856

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LL37 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Banded Krait

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. Use at a dilution of 1:100 to 1: 2000. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide from Bungarus fasciatus Cathelicidin antimicrobial peptide (CAP-18, antibacterial protein LL-37, CAMP, CRAMP) conjugated to an immunogenic carrier protein has been used as the antigen.
Immunogen Region	
Specificity	Specific for Cathelicidin antimicrobial peptide.
Immunogen Sequence	