

Anti-Nicotinic acetylcholine R alpha-10 antibody (C-Term) (STJ13100905)

STJ13100905

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nicotinic acetylcholine R alpha-10 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	A synthetic peptide from c-terminal region of human Neuronal acetylcholine receptor subunit alpha-10 conjugated to blue carrier protein was used as the antigen. The antigen shares 85% identity with mouse and rat sequences.
Immunogen Region	C-Term
Specificity	Specific for CHRNA10.
Immunogen Sequence	