

Anti-CNTNAP2 antibody (ECD) (STJ13100906)

STJ13100906

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CNTNAP2 (ECD) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	26047
Gene Symbol	CNTNAP2
Uniprot ID	CNTP2_HUMAN
Immunogen	A synthetic peptide from extracellular domain of human Contactin associated protein like 2 (CNTNAP2) conjugated to blue carrier protein was used as the antigen.
Immunogen Region	ECD
Specificity	Specific for CNTNAP2.
Immunogen Sequence	