

Anti-Contactin 2 antibody (650-700 aa) (STJ13100991)

STJ13100991

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Contactin 2 (650-700 aa) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|--|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | IgG purified |
| Dilution Range | IHC, WB (confirmed by recombinant protein). A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications. |
| Formulation | Lyophilised |
| Isotype | IgG |
| Storage | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. |
| Instruction | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | A synthetic peptide from aa region 650-700 of human Contactin2 conjugated to an immunogenic carrier protein was used as the antigen. |
| Immunogen Region | 650-700 aa |
| Specificity | Specific for contactin-2. |
| Immunogen Sequence | |