

Anti-EXOC4 antibody (N-Term) (STJ13101070)

STJ13101070

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-EXOC4 (N-Term) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	20336
Gene Symbol	Exoc4
Uniprot ID	EXOC4_MOUSE
Immunogen	A synthetic peptide from the n-terminal of human Exocyst complex component 4 (EXOC4, SEC8, SEC8 L1) conjugated to blue carrier protein was used as the antigen. The antigen is homologous in mouse and rat.
Immunogen Region	N-Term
Specificity	Specific for EXOC4.
Immunogen Sequence	