

Anti-FGF13 antibody (2-50 aa) (STJ13101111)

STJ13101111

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-FGF13 (2-50 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	2258
Gene Symbol	FGF13
Uniprot ID	FGF13_HUMAN
Immunogen	A synthetic peptide from aa region 2-50 of mouse FGF13 conjugated to an immunogenic carrier protein was used as the antigen. The peptide is homologous in human and rat.
Immunogen Region	2-50 aa
Specificity	Specific for FGF13.
Immunogen Sequence	