

Anti-GNB3 antibody (2-50 aa) (STJ13101132)

STJ13101132

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GNB3 (2-50 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 1000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	2784
Gene Symbol	GNB3
Uniprot ID	GBB3_HUMAN
Immunogen	A synthetic peptide from aa region 2-50 of human GNB3 conjugated to blue carrier protein was used as the antigen. The peptide is homologous in rat and mouse.
Immunogen Region	2-50 aa
Specificity	Specific for GBB3.
Immunogen Sequence	