

Anti-Phospho-LIMK2-S291/S293 antibody (STJ13101311)

STJ13101311

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-Phospho-LIMK2-S291/S293 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. This antibody may also recognise the unphosphorylated form. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	16886
Gene Symbol	Limk2
Uniprot ID	LIMK2_MOUSE
Immunogen	A synthetic peptide from the phospho S291 & S293 region of human LIMK2 (LIM domain kinase 2, LIMK-2) conjugated to an immunogenic carrier protein was used as the immunogen. The antigen is homologous in rat, mouse.
Immunogen Region	
Specificity	Specific for LIMK2.
Immunogen Sequence	