

Anti-LATS2 antibody (5-50 aa) (STJ13101313)

STJ13101313

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-LATS2 (5-50 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	50523
Gene Symbol	Lats2
Uniprot ID	LATS2_MOUSE
Immunogen	A synthetic peptide from aa region 5-50 of mouse LATS2 conjugated to blue carrier protein was used as the antigen.
Immunogen Region	5-50 aa
Specificity	Specific for LATS2.
Immunogen Sequence	