

Anti-Lamp1 antibody (STJ13101317)

STJ13101317

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lamp1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	16783
Gene Symbol	Lamp1
Uniprot ID	LAMP1_MOUSE
Immunogen	A synthetic peptide to from second luminal domain of mouse Lamp1 (CD107a) conjugated to blue carrier protein was used as the antigen. The peptide is homologous in rat.
Immunogen Region	
Specificity	Specific for Lamp1.
Immunogen Sequence	