

Anti-p75NTR antibody [ME20.4] (STJ13101332)

STJ13101332

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-p75NTR is suitable for use in Immunohistochemistry, Western Blot and Flow Cytometry research applications.
Applications	IHC/WB/FC
Host/Source	Mouse
Reactivity	Human/Marmoset/Dog/Cat/Pig/Rabbit

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ME20.4
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB, FC. A concentration of 2-10 $\mu\text{g/ml}$ is recommended. The optimal working concentration should be determined by the end user. This antibody may also be used for immunoabsorptive purification of NGF-receptors. See also the eLab Book at the end o
Formulation	Lyophilised
Isotype	IgG1
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Human melanoma cell lines
Immunogen Region	
Specificity	Specific for p75NTR.
Immunogen Sequence	