

Anti-MAP1LC3B antibody (STJ13101341)

STJ13101341

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-MAP1LC3B is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC, WB. Use at a dilution of 1:200 to 1:2000. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	81631
Gene Symbol	MAP1LC3B
Uniprot ID	MLP3B_HUMAN
Immunogen	A synthetic peptide from human LC3 B conjugated to blue carrier protein has been used as the antigen.
Immunogen Region	
Specificity	Specific for MAP1LC3B.
Immunogen Sequence	