

Anti-Mucolipin 2 antibody (C-Term) (STJ13101344)

STJ13101344

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Sheep polyclonal antibody anti-Mucolipin 2 (C-Term) is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Sheep |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | Affinity purification |
| Dilution Range | IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. |
| Formulation | Lyophilised |
| Isotype | |
| Storage | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. |
| Instruction | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | |
| Gene Symbol | |
| Uniprot ID | |
| Immunogen | A synthetic peptide from the c-terminal region of human Mucolipin 2 conjugated to an immunogenic carrier protein was used as the antigen. |
| Immunogen Region | C-Term |
| Specificity | Specific for Mucolipin 2. |
| Immunogen Sequence | |