

Anti-NAV3 antibody (Internal Region) (STJ13101475)

STJ13101475

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NAV3 (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	89795
Gene Symbol	NAV3
Uniprot ID	NAV3_HUMAN
Immunogen	A synthetic peptide from internal region of human Neuron navigator 3 (NAV3, Steerin-3) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	Internal Region
Specificity	Specific for NAV3.
Immunogen Sequence	