

Anti-NBR1 antibody (290-340 aa) (STJ13101506)

STJ13101506

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-NBR1 (290-340 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse/Rat/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	17966
Gene Symbol	Nbr1
Uniprot ID	NBR1_MOUSE
Immunogen	A synthetic peptide from aa region 290-340 of mouse NBR1 conjugated to an immunogenic carrier protein was used as the antigen. The peptide is homologous in rat and shares 98% identity with human sequence.
Immunogen Region	290-340 aa
Specificity	Specific for NBR1
Immunogen Sequence	