

Anti-OR7D4 antibody (STJ13101549)

STJ13101549

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-OR7D4 is suitable for use in Immunohistochemistry and Western Blot research applications. |
| Applications | IHC/WB |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|-----------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | IgG purified |
| Dilution Range | IHC, WB. A dilution of 1:200 to 1:2000 is recommended. The optimal dilution should be determined by the end user. |
| Formulation | Lyophilised |
| Isotype | IgG |
| Storage | Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term. |
| Instruction | When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 125958 |
| Gene Symbol | OR7D4 |
| Uniprot ID | OR7D4_HUMAN |
| Immunogen | A synthetic peptide from human olfactory receptor (OR7D4, OR19-7, OR19-B) conjugated to blue carrier protein was used as the antigen. |
| Immunogen Region | |
| Specificity | Specific for Olfactory receptor OR7D4. |
| Immunogen Sequence | |