

Anti-KAT2B antibody (400-460 aa) (STJ13101687)

STJ13101687

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-KAT2B (400-460 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	18519
Gene Symbol	Kat2b
Uniprot ID	KAT2B_MOUSE
Immunogen	A synthetic peptide from aa region 400-460 of mouse PCAF conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	400-460 aa
Specificity	Specific for PCAF
Immunogen Sequence	