

Anti-TAAR1 antibody (210-260 aa) (STJ13101748)

STJ13101748

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAAR1 (210-260 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1: 300 to 1: 2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	134864
Gene Symbol	TAAR1
Uniprot ID	TAAR1_HUMAN
Immunogen	A synthetic peptide from aa region 210-260 of human TAAR1 conjugated to blue carrier Protein was used as the antigen.
Immunogen Region	210-260 aa
Specificity	Specific for TAAR1.
Immunogen Sequence	