

Anti-MRGPRG antibody (STJ13101769)

STJ13101769

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-MRGPRG is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	386746
Gene Symbol	MRGPRG
Uniprot ID	MRGRG_HUMAN
Immunogen	A synthetic peptide from the 4th cytoplasmic loop of human Mas-related G-protein coupled receptor member G (MRGPRG) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for MRGPRG.
Immunogen Sequence	