

Anti-TAAR5 antibody (C-Term) (STJ13101778)

STJ13101778

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TAAR5 (C-Term) is suitable for use in Immunohistochemistry research applications.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	215854
Gene Symbol	Taar5
Uniprot ID	TAAR5_MOUSE
Immunogen	A synthetic peptide from c-terminal region of mouse trace amine-associated receptor 5 (TAAR5) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	C-Term
Specificity	Specific for TAAR5.
Immunogen Sequence	