

Anti-TRPM2 antibody (4th Cytoplasmic Domain) (STJ13101783)

STJ13101783

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TRPM2 (4th Cytoplasmic Domain) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	7226
Gene Symbol	TRPM2
Uniprot ID	TRPM2_HUMAN
Immunogen	A synthetic peptide from the 4th cytoplasmic domain of human TRPM2 (TRPC7, LTRPC2) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	4th Cytoplasmic Domain
Specificity	Specific for TRPM2.
Immunogen Sequence	