

Anti-TMC2 antibody (STJ13101848)

STJ13101848

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-TMC2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB (confirmed by recombinant protein). A dilution of 1:200 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	192140
Gene Symbol	Tmc2
Uniprot ID	TMC2_MOUSE
Immunogen	A synthetic peptide from the 1 st cytoplasmic loop of mouse transmembrane channel-like protein 2 (TMC2) conjugated to an immunogenic carrier protein was used as the immunogen.
Immunogen Region	
Specificity	Specific for TMC2.
Immunogen Sequence	