

Anti-OR2AG2 antibody (STJ13101880)

STJ13101880

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OR2AG2 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:200 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	338755
Gene Symbol	OR2AG2
Uniprot ID	O2AG2_HUMAN
Immunogen	A synthetic peptide from human olfactory receptor (OR2AG2) conjugated to an immunogenic carrier protein was used as the antigen.
Immunogen Region	
Specificity	Specific for olfactory receptor OR2AG2.
Immunogen Sequence	