

**Anti-OR51E1 antibody (235-275 aa) (STJ13101893)**

STJ13101893

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Sheep polyclonal antibody anti-OR51E1 (235-275 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC/WB
<b>Host/Source</b>	Sheep
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user. Not yet tested in other applications.
<b>Formulation</b>	Lyophilised
<b>Isotype</b>	
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	143503
<b>Gene Symbol</b>	OR51E1
<b>Uniprot ID</b>	O51E1_HUMAN
<b>Immunogen</b>	A synthetic peptide from aa region 235-275 of human OR51E1 conjugated to blue carrier protein was used as the antigen.
<b>Immunogen Region</b>	235-275 aa
<b>Specificity</b>	Specific for OR51E1.
<b>Immunogen Sequence</b>	