

Anti-OR2Y1 antibody (STJ13101896)

STJ13101896

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OR2Y1 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:3000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	134083
Gene Symbol	OR2Y1
Uniprot ID	OR2Y1_HUMAN
Immunogen	A synthetic peptide as a part of human olfactory receptor OR2Y1 (OR5-2) conjugated to an immunogenic carrier protein was used as the immunogen.
Immunogen Region	
Specificity	Specific for Olfactory receptor OR2Y1.
Immunogen Sequence	