

Anti-OR3A3 antibody (STJ13101897)

STJ13101897

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OR3A3 is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	8392
Gene Symbol	OR3A3
Uniprot ID	OR3A3_HUMAN
Immunogen	A synthetic peptide from human olfactory receptor OR3A3 (OR3A7, OR3A8P, OR3A6, Olfactory receptor 17-201, OR17-201, Olfactory receptor 3A3, OR17-16, OR17-137, OR17-201) conjugated to an immunogenic carrier protein was used as the antigen. The peptide
Immunogen Region	
Specificity	Specific for Olfactory receptor OR3A3.
Immunogen Sequence	