

Anti-OR4C11 antibody (220-260 aa) (STJ13101900)

STJ13101900

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-OR4C11 (220-260 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A dilution of 1:300 to 1:2000 is recommended. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	219429
Gene Symbol	OR4C11
Uniprot ID	OR4CB_HUMAN
Immunogen	A synthetic peptide from aa region 220-260 of human olfactory receptor (OR4C11) conjugated to blue carrier protein was used as the antigen.
Immunogen Region	220-260 aa
Specificity	Specific for Olfactory receptor OR4C11.
Immunogen Sequence	