

Anti-RhoB antibody (STJ13101905)

STJ13101905

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RhoB is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	IgG purified
Dilution Range	IHC, WB. A concentration of 10-50 µg/ml is recommended. The optimal concentration should be determined by the end user.
Formulation	Lyophilised
Isotype	IgG
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	388
Gene Symbol	RHOB
Uniprot ID	RHOB_HUMAN
Immunogen	A synthetic peptide as a part of human, rat, mouse Rho-related GTP-binding protein RhoB (ARH6, ARHB, H6) conjugated to an immunogenic carrier protein was used as the immunogen.
Immunogen Region	
Specificity	Specific for RhOB.
Immunogen Sequence	