

## Anti-RAB5A antibody (Internal Region) (STJ13101919)

STJ13101919

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-RAB5A (Internal Region) is suitable for use in Immunohistochemistry and Western Blot research applications.
<b>Applications</b>	IHC/WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	IgG purified
<b>Dilution Range</b>	IHC, WB. Use at a concentration of 10-50 µg/ml. The optimal dilution should be determined by the end user. Not yet tested in other applications.
<b>Formulation</b>	Lyophilised
<b>Isotype</b>	IgG
<b>Storage</b>	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
<b>Instruction</b>	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">271457</a>
<b>Gene Symbol</b>	<a href="#">Rab5a</a>
<b>Uniprot ID</b>	<a href="#">RAB5A_MOUSE</a>
<b>Immunogen</b>	A synthetic peptide from the internal region of mouse Ras-related protein Rab-5A (RAB5A) conjugated to an immunogenic carrier protein has been used as the immunogen. The antigen is homologous in rat.
<b>Immunogen Region</b>	Internal Region
<b>Specificity</b>	Specific for RAB5A.
<b>Immunogen Sequence</b>	