

Anti-SQSTM1 antibody (300-360 aa) (STJ13101971)

STJ13101971

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Sheep polyclonal antibody anti-SQSTM1 (300-360 aa) is suitable for use in Immunohistochemistry and Western Blot research applications.
Applications	IHC/WB
Host/Source	Sheep
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC (both on frozen and paraffin embedded tissues) , WB. Use at a dilution of 1:200 to 1:2000. The optimal dilution should be determined by the end user.
Formulation	Lyophilised
Isotype	
Storage	Maintain the lyophilised/reconstituted antibodies frozen at -20°C for long term storage and refrigerated at 2-8°C for a shorter term.
Instruction	When reconstituting, Glycerol (1:1) may be added for an additional stability. Avoid freeze and thaw cycles.

TARGET INFORMATION

Gene ID	8878
Gene Symbol	SQSTM1
Uniprot ID	SQSTM_HUMAN
Immunogen	A synthetic peptide from aa region 300-360 of human SQSTM1 has been used as the antigen.
Immunogen Region	300-360 aa
Specificity	Specific for sequestosome-1.
Immunogen Sequence	